

“Despite improvements in some regions, water pollution is on the rise globally. And unless substantial progress is made in regulation and enforcement, pollution is expected to increase as a result of economic development driven by urbanization, industries and intensive agriculture systems.”

Source: 3rd United Nations World Water Development Report: WATER IN A CHANGING WORLD



## WATER ANALYSIS

Water is among the most precious natural resources. Every living thing on earth needs water to survive and it is of vitally importance to every aspect of our lives. Humanity has always been aware of the importance of water. Aeschylus, more than two thousand years ago, said: “By polluting clear water with slime, you will never find good drinking water”.

But nowadays each water source can contain dramatically different levels of pollution. Measurements of pollutant levels in water provide the most fundamental indicator of the status of this resource and are critical and meaningful evaluation of the quality of water.

Monitoring of contaminants is a matter of great importance as it influences human and environmental health.

In the following pages a number of solutions that allow achieving accurate, precise, and reliable results to match the requirements of the regulated water analysis are presented.

DANI has continued to implement innovative techniques by introducing in the market very reliable, extremely versatile and easy to use instruments that meet water monitoring real requirements.

Decades of experience in developing applications for the Environmental industry are the basis of the Master DWA DANI Water Analyzers.

# WATER ANALYSIS



*Time is Money.*



Benjamin Franklin

Time is always a key factor in today's laboratories productivity.

## Master your Time with the Master DANI Water Analyzers.

The ability to provide the proper configuration to meet the most challenging analytical demands comes from a long and proven experience and a deep industry knowledge. As requirements are constantly changing, even a highly reliable instrumentation could not be enough to succeed in getting trustworthy results: complete and guaranteed solutions are essential to comply with the latest industry standards and specifications.

After a long working relationship with its customers to know and to best match their real needs, DANI Instruments has developed key analytical solutions that cover a broad array of applications, requirements and protocols in the environmental industry.

Master DANI Water Analyzers are **PRE-CONFIGURED, PRE-ASSEMBLED AND FACTORY-TESTED SYSTEMS** specifically designed for specific analyses. The analyzers include the **HARDWARE**, the **SOFTWARE**, **COLUMNS AND CONSUMABLES**, the **OPTIMIZED ANALYSIS METHOD**, the **ANALYTICAL CONDITIONS**, and the **DOCUMENTATION** to run up your analysis from day one.

### PRE-CONFIGURED, PRE-ASSEMBLED AND FACTORY-TESTED SYSTEMS

The installation process is faster than ever before and all the startup procedure is oversimplified ensuring immediate analytical performance and results.

### HARDWARE AND SOFTWARE

DWA Analyzers are pre-engineered systems based on the versatility, flexibility and robustness of the proven Master GC hardware. All the Master GC parameters are set prior the shipment.

### COLUMNS AND CONSUMABLES

No more doubts about the proper column, parts and supplies. DWA Analyzers are delivered with all you may need for your analysis.\*

### OPTIMIZED ANALYSIS METHOD

Analytical methods are pre-loaded to be immediately used for the determination of pollutants in water. Whenever possible, reduction of analysis time and amounts of toxic solvents are considered. Method development time and costs are thus dramatically reduced.

### ANALYTICAL CONDITIONS

DWA Analyzers are designed to perfectly accomplish the analytical conditions of interest.

### DOCUMENTATION

A getting started manual, calibration and method files, and all the information for a quick startup are included.



### View your Analyzer:

**Master DWA-114**  
Halogenated Hydrocarbons in Drinking Water Analyzer

**Master DWA-115**  
Volatile Aromatic Compounds in Drinking Water Analyzer

**Master DWA-119**  
Volatile Aromatic Hydrocarbons in Water Analyzer

**Master DWA-120**  
VOCs in Water Analyzer

**Master DWA-121**  
PCBs in Water Analyzer

**Master DWA-049**  
BTEX and Styrene in Water Analyzer

**Master DWA-122**  
Fuel Oxygenates in Water Analyzer

**Master DWA-123**  
Hydrocarbon Oil Index in Water Analyzer

\* Chemicals are not supplied



**Master GC**  
**Fast Gas Chromatograph**  
*High Productivity*  
*Accuracy and Precision*  
*Flexibility and Upgradeability*  
*User Friendly Interface*



The versatile and flexible Master GC delivers unsurpassed analytical capabilities meeting today's laboratories productivity requirements. The GC was uniquely designed to perform conventional and fast gas chromatographic analyses. The primary goal of Fast GC is to maintain proper resolving power in shorter analysis run times by using adequate instrumentation and analytical columns in combination with optimized method parameters. The Master GC features a maximum heating rate up to 140° C/min and a typical cooling time of 4 min. In addition, DANI offers a variety of detectors engineered with fast electronics to handle sharp peaks; data acquisition rates of up to 300 Hz are performed.


The Master TOF-MS detector performs the fastest acquisition rates (1000 spectra/s) and the widest linear dynamic range ( $10^5$ ) available on the market. These capabilities are offered in an extremely compact bench-top instrument. In combination with the Master GC, the system is the ideal solution for Fast GC and GCxGC laboratories.

The Master LAB Software offers the proper tool for the reliable control of the system, from autotuning procedures to GC and sample sequence management. An original deconvolution algorithm capable to handle a large amount of information in a smart and effective way provides trustworthy identification of trace compounds even in complex matrices.



**Master TOF**  
**Time of Flight GC/MS**  
*Extremely Compact Design*  
*High Productivity*  
*Powerful Software Solution*  
*Walkaway Automation*

**Master SHS**  
**Static Headspace Sampler**  
*A Robust and Flexible System*  
*to Meet Complex and Versatile*  
*Needs*



The Master SHS delivers the highest performances to overcome daily new challenges and supplies trustworthy and enhanced results. The highest sample capacity and the unlimited priority sample position provide straightforward results for virtually any analytical need in real time. The Valve&Loop Technique, the known and fixed volume of the sample, the accurate temperature control and the entirely chemically inert sample flow path guarantee outstanding repeatability and avoid the risk of false results, sample loss or recondensation.

The Master DHS/P&T provides the most versatile, state-of-the-art system for headspace analysis featuring the capabilities of a Purge&Trap system. It combines the high sensitivity of the Dynamic Headspace technique with the productivity, ease of use, and flexibility of a completely automated solution. The Master DHS/P&T offers up to a 100-fold increase in sensitivity over conventional headspace techniques and assures detection limits beyond capability of SPME.



**Master DHS/P&T**  
**Dynamic Headspace and**  
**Purge&Trap Sampler**  
*A Dynamic Approach*  
*to High Sensitivity*  
*Headspace Analysis*

**Master TD**  
**Thermal Desorber**  
*The Ultimate Solution for*  
*High Sensitivity Detection of*  
*Volatiles*



The Master TD offers superior sensitivity, versatility, and productivity for the extraction of volatile and semi-volatile compounds from air and solid matrices. The excellent analytical performances of the system is guaranteed by the two-stage thermal desorption process and supported by the patented "Instant Desorption" of the trap. This design assures the complete transfer of the analytes and their injection into the analytical column in a narrow band to preserve chromatographic resolution and accuracy. The fully automated control of the system provides high sample capacity and optimal sampling tube processing for maximum system productivity.

# WATER ANALYSIS

## TABLE OF CONTENTS

### Halogenated Hydrocarbons

1,1-dichloroethylene  
methylene chloride  
*trans*-1,2-dichloroethylene  
1,1-dichloroethane  
2,2-dichloropropane  
*cis*-1,2-dichloroethylene  
chloroform  
bromochloromethane  
1,1,1-trichloroethane  
1,1-dichloropropylene  
carbon tetrachloride  
1,2-dichloroethane  
trichloroethylene  
1,2-dichloropropane  
bromodichloromethane  
dibromomethane  
*cis*-1,3-dichloropropylene  
*trans*-1,3-dichloropropylene  
1,1,2-trichloroethane

1,3-dichloropropane  
tetrachloroethylene  
dibromochloromethane  
1,2-dibromoethane  
chlorobenzene  
1,1,1,2-tetrachloroethane  
bromoform  
1,1,2,2-tetrachloroethane  
1,2,3-trichloropropane  
bromobenzene  
2-chlorotoluene  
4-chlorotoluene  
1,3-dichlorobenzene  
1,4-dichlorobenzene  
1,2-dichlorobenzene  
1,2-dibromo-3-chloropropane  
1,2,4-trichlorobenzene  
1,2,3-trichlorobenzene

EPA Method 502.2

Application Note AN 114  
Halogenated Hydrocarbons in Drinking Water

*Master* DWA-114  
Halogenated Hydrocarbons in Water Analyzer

### Volatile Aromatic Organic Compounds

benzene  
toluene  
ethylbenzene  
*m*-xylene  
*p*-xylene  
*o*-xylene  
styrene  
iso-propylbenzene  
*n*-propylbenzene  
bromobenzene  
1,3,5-trimethylbenzene

2-chlorotoluene  
4-chlorotoluene  
*tert*-butylbenzene  
1,2,4-trimethylbenzylene  
*sec*-Butylbenzene  
*p*-Isopropyltoluene  
1,3-Dichlorobenzene  
1,4-Dichlorobenzene  
*n*-butylbenzene  
1,2-Dichlorobenzene

EPA Method 502.2

Application Note AN 115  
Volatile Aromatic Organic Compounds in Water

*Master* DWA-115  
Volatile Aromatic Organic Compounds in Water Analyzer

### Volatile Aromatic Hydrocarbons

benzene  
toluene  
ethylbenzene  
*p*-Xylene  
*m*-Xylene  
*o*-Xylene

Application Note AN 119  
Dynamic Determination of Volatile Hydrocarbons

*Master* DWA-119  
Volatile Hydrocarbons in Water Analyzer

### VOCs

1,1-dichloroethene  
methylene chloride  
*trans*-1,2-dichloroethene  
*cis*-1,2-Dichloroethene  
chloroform  
carbon tetrachloride  
1,1,1-trichloroethane  
benzene  
1,2-dichloroethane  
trichloroethene  
1,2-dichloropropane  
bromodichloromethane

*cis*-1,3-dichloropropene  
toluene  
*trans*-1,3-dichloropropene  
tetrachloroethane  
1,1,2-trichloroethane  
dibromochloromethane  
*m*-Xylene  
*p*-Xylene  
*o*-Xylene  
bromoform  
1,4-dichlorobenzene

Application Note AN 120  
Dynamic and Purge & Trap Determination of VOCs

*Master* DWA-120  
VOCs in Water Analyzer

BTEX and Styrene

benzene  
toluene  
ethylbenzene  
*p*-xylene  
*m*-xylene  
*o*-xylene  
styrene

Application Note AN 049  
Determination of BTEX and Styrene in Water

Master DWA-049  
BTEX and Styrene in Water Analyzer

Phenols

phenol  
2-chlorophenol  
2,4-dimethylphenol  
2,4-dichlorophenol  
4-chloro-3-methylphenol

2,4,6-trichlorophenol  
2,4-dinitrophenol  
4-nitrophenol  
2-methyl-4-dinitrophenol  
pentachlorophenol

EPA Method 604  
EPA Method 8041

Application Note AN 007  
Phenols: Fast GC Analysis

Poly Aromatic Hydrocarbons (PAHs)

naphtalene  
2-methylnaphtalene  
1-methylnaphtalene  
acenaphthylene  
acenaphtene  
fluorene  
phenanthrene  
anthracene  
fluoranthene

pyrene  
benzo[a]anthracene  
chrysene  
benzo[b]fluoranthene  
benzo[k]fluoranthene  
benzo[a]pyrene  
indeno[1,2,3-cd]pyrene  
dibenzo[a,h]anthracene  
benzo[g,h,i]perylene

EPA Method 8100

Application Note AN 003  
Poly Aromatic Hydrocarbons (PAHs) - Fast GC Application

PCBs

2,4,5,6-tetrachloro-*m*-xylene (S.S)  
2,3-dichlorobiphenyl  
2,2',5'-trichlorobiphenyl  
2,4',5'-trichlorobiphenyl  
2,2',5,5'-tetrachlorobiphenyl  
2,2',3,5'-tetrachlorobiphenyl  
2,3',4,4'-tetrachlorobiphenyl  
2,2',4,5,5'-pentachlorobiphenyl  
2,2',3,4,5'-pentachlorobiphenyl  
2,3,3',4',6'-pentachlorobiphenyl

2,2',3,5,5',6-hexachlorobiphenyl  
2,2',4,4',5,5'-hexachlorobiphenyl  
2,2',3,4,5,5'-hexachlorobiphenyl  
2,2',3,4,4',5'-hexachlorobiphenyl  
2,2',3,4',5,5',6-heptachlorobiphenyl  
2,2',3,4,4',5',6-heptachlorobiphenyl  
2,2',3,4,4',5,5'-heptachlorobiphenyl  
2,2',3,3',4,4',5'-heptachlorobiphenyl  
2,2',3,3',4,4',5,5'-nonachlorobiphenyl  
decachlorobiphenyl (I.S.)

EPA Method 8082

Application Note AN 121  
Fast GC Approach for PCBs Determination

Master DWA-121  
PCBs in Water Analyzer

Fuel Oxygenates

TBA  
MtBE  
DIPE  
EtBE  
TAME

Application Note AN 122  
Fuel Oxygenates in Water

Master DWA-122  
Fuel Oxygenates in Water Analyzer

Hydrocarbon Oil Index C7-C40

C7C16C25C34  
C8C17C26C35  
C9C18C27C36  
C10C19C28C37  
C11C20C29C38  
C12C21C30C39  
C13C22C31C40  
C14C23C32  
C15C24C33

ISO 9377-2 Method

Application Note AN 123  
Hydrocarbon Oil Index in Water

Master DWA-123  
Hydrocarbon Oil Index in Water Analyzer

# WATER ANALYSIS

## TABLE OF CONTENTS

### Phthalates

dimethyl-phtalate  
diethyl-phtalate  
di-*n*-butyl phtalate

butyl benzyl phtalate  
*bis* (2-ethylhexyl) phtalate  
di-*n*-octyl phtalate

Application Note AN 052  
Determination of Phthalates

### Haloethers

*bis*(2-chloroethyl) ether  
*bis*(2-chloroisopropyl)ether  
*bis*(2-chloroethoxy)ether  
4-chlorophenylphenyl ether  
4-bromophenyl phenyl ether

EPA Method 611/8111

Application Note AN 069  
Determination of Haloethers

### Nitrosamines

*n*-nitrosodimethylamine  
*n*-nitrosodi-*n*-propylamine  
*n*-nitrosodiphenylamine

Application Note AN 067  
Determination of Nitrosamines

### Phenols and Chlorophenols

phenol  
2-chlorophenol  
2-nitrophenol  
2,4-dichlorophenol  
4-chloro-3-methylphenol

2,4,6-trichlorophenol  
2,4-dinitrophenol  
4-nitrophenol  
2-methyl-4,6-ditrophenol  
pentachlorophenol

EPA Method 604

Application Note AN 066  
Determination of Phenols and Chlorophenols

### Organochlorinated Pesticides

$\alpha$  - BHC  
 $\beta$  - BHC  
 $\gamma$  - BHC  
 $\delta$  - BHC  
heptachlor  
aldrin  
heptachlor epoxide  
endosulfan I

4,4' DDE  
dieldrin  
endrin  
4,4' DDD  
endosulfan II  
endrin aldehyde  
4,4' DDT  
endosulfan sulfate

EPA Method 608/8081

Application Note AN 063  
Determination of Organochlorinated Pesticides

### Organophosphorus Pesticides

trichlorfon  
tionazine  
phorate  
diazinone  
CH<sub>3</sub>-parathion  
chlorpyrifos  
parathion  
quinalphos  
metidathion  
yrithion  
C<sub>2</sub>H<sub>5</sub>-azinphos

phosdrin  
ethoprophos  
phonophos  
CH<sub>3</sub>-chlorpyrifos  
CH<sub>3</sub>-pirimiphos  
malathion  
pirimiphos  
C<sub>2</sub>H<sub>5</sub>-bromophos  
ethion  
CH<sub>3</sub>-azinphos

Application Note AN 095  
Determination of Organophosphorus Pesticides





# WATER ANALYSIS

## Application Notes

# WATER ANALYSIS

## Halogenated Hydrocarbons in Drinking Water

### EPA Method 502.2

#### Application Note AN 114

Halogenated hydrocarbons can be found in appreciable amounts in surface and drinking waters. Often these contaminants are the result of the chlorination of raw, groundwater or wastewater in order to achieve drinking water quality standards. Contamination of water by halogenated hydrocarbons has been disclosed to be toxic for humans. It's therefore of primary importance to have an easy-to-use, **RELIABLE** and **COMPLETELY AUTOMATED** method for this type of analysis. In the following application Gas Chromatography with Purge&Trap concentration and Electron Capture Detection has been applied to the determination of halogenated hydrocarbons in water, **ACCORDING TO EPA METHOD 502.2 REQUIREMENTS**.

*Do you wish to streamline the process, boost your productivity and save time and hassles?*

DANI Water Analyzer DWA-114 is the ready-to-go solution to attain the maximum performance in the shortest time for your analysis of Halogenated Hydrocarbons in Drinking Water.

**Master DWA-114**

#### COMPLETE AUTOMATION OF ALL THE ANALYTICAL STEPS FOR AN INCREASED PRODUCTIVITY

Master DHS/P&T allows sample overlapping: the system automatically controls that the next sample is thermostatted during the GC analysis of the previous one. Solutions with different concentrations can also be prepared in a fully automatic way.

#### RELIABLE SYSTEM WITH NO CARRY-OVER RISK

Each sample is placed in a disposable 20-mL headspace vial. No additional workload of cleaning glassware or line purging is necessary, and the sample needle undergoes on automated cleaning cycle during the baking phase. Furthermore, the entire sample flow path is chemically inert. All these features prevent carry-over effects.



### Analysis Conditions

#### Master GC Parameters:


Oven	35°C (8 min), 4°C/min, 240°C (1 min)
Detector	ECD 300°C - 40ml/min N <sub>2</sub>
Injector	SL/IN (220°C)
Carrier	Helium, 3.5 mL/min, (split 1:2)
Column	Vocol 60 m x 0.32 mm i.d. x 3µm d <sub>f</sub>

#### Master DHS/P&T Parameters - Purging Mode

Incubation	60°C
Stripping	3 min, 120 mL/min, Trap -10°C
Injection	3 min, Dew Stop 0°C, Trap 295°C
Baking	10 min, 80 mL/min Trap 300°C, Dew Stop 200°C
Transfer Line	250°C
Switching Valve	250°C
Trap Material	Tenax/Carbotrap/Carbosieve
Sample Volume	10mL



## Halogenated Hydrocarbons in Drinking Water EPA Method 502.2



	COMPOUNDS	RSD%	MDL	R	RANGE	EPA 502.2 RSD%	EPA 502.2 MDL
			ppb		ppb		ppb
1	1,1-dichloroethylene	1.5	0.0100	0.9570	0.1000-2	2.8	0.07
2	methylene chloride	2.0	0.0020	0.9960	0.0100-10	2.9	0.02
3	<i>trans</i> -1,2-dichloroethylene	3.0	0.0300	0.9970	0.0100-10	3.7	0.06
4	1,1-dichloroethane	2.2	0.0030	0.9970	0.0200-2	5.7	0.07
5	2,2-dichloropropane	1.1	0.0050	0.9990	0.0002-10	3.4	0.05
6	<i>cis</i> -1,2-dichloroethylene	16.6	0.0100	0.9990	0.0200-10	3.3	0.01
7	chloroform	0.8	0.0005	0.9850	0.0400-2	2.5	0.02
8	bromochloromethane	0.8	0.0005	0.9940	0.0400-1	3.0	0.01
9	1,1,1-trichloroethane	0.9	0.0005	0.9980	0.0020-2	3.3	0.03
10	1,1-dichloropropylene	0.7	0.0003	0.9940	0.0100-1	3.3	0.02
10	carbon tetrachloride	0.7	0.0003	0.9970	0.0020-1	3.6	0.01
11	1,2-dichloroethane	0.5	0.0020	0.9980	0.0100-1	3.8	0.03
12	trichloroethylene	1.2	0.0005	0.9980	0.0100-1	3.6	0.01
13	1,2-dichloropropane	1.3	0.0020	0.9960	0.0020-1	3.7	0.01
14	dromodichloromethane	0.8	0.0003	0.9950	0.040-0.1	2.9	0.02
15	dibromomethane	1.6	0.0003	0.9640	0.0100-1	1.5	0.02
16	<i>cis</i> -1,3-dichloropropylene	2.0	0.0007	0.9940	0.0002-2	3.7	0.06
17	<i>trans</i> -1,3-dichloropropylene	1.6	0.0010	0.9970	0.0002-2	33.7	0.01
18	1,1,2-trichloroethane	1.5	0.0014	0.9950	0.0002-2	5.6	N.D.
19	1,3-dichloropropane	0.7	0.0004	0.9980	0.0400-2	3.1	0.03
19	tetrachloroethylene	0.7	0.0004	0.9980	0.0400-2	2.5	0.04
20	dibromochloromethane	0.5	0.0003	0.9600	0.0400-1	2.8	0.08
21	1,2-dibromoethane	0.9	0.0005	0.9940	0.0020-1	6.7	2.20
22	chlorobenzene	0.9	0.0003	0.9950	0.0100-1	3.6	0.01
22	1,1,1,2-tetrachloroethane	0.9	0.0003	0.9950	0.0100-1	2.3	0.01
23	bromoform	1.0	0.0005	0.9980	0.0100-1	5.2	1.60
24	1,1,2,2-tetrachloroethane	1.6	0.0007	0.9987	0.0002-1	6.8	0.01
25	1,2,3-trichloropropane	0.8	0.0025	0.9972	0.0100-1	2.3	0.40
26	bromobenzene	1.6	0.0120	0.9951	0.2000-2	2.7	0.03
27	2-chlorotoluene	3.3	0.0020	0.9977	0.0200-2	2.7	0.01
27	4-chlorotoluene	3.3	0.0020	0.9976	0.0200-2	3.2	0.01
28	1,3-dichlorobenzene	0.2	0.0070	0.9993	0.2000-10	4.0	0.02
29	1,4-dichlorobenzene	3.4	0.0180	0.9989	0.4000-10	2.3	0.01
30	1,2-dichlorobenzene	3.0	0.0100	0.9985	0.2000-10	1.5	0.02
31	1,2-dibromo-3-chloropropane	4.0	0.0050	0.9976	0.2000-2	11.3	3.00
32	1,2,4-trichlorobenzene	2.8	0.0040	0.9989	0.1000-10	2.1	0.03
33	1,2,3-trichlorobenzene	2.5	0.0011	0.9994	0.1000-10	3.1	0.03

# WATER ANALYSIS

## Volatile Aromatic Organic Compounds in Drinking Water A Dynamic Approach to EPA Method 502.2

### Application Note AN 115

Drinking water containing high level of volatile organic compounds may be harmful to human health. The U.S. Environmental Protection Agency estimates that VOCs are present in one fifth of the Nation's water supplies. Some VOCs are mutagens, teratogens, and carcinogens.

EPA Method 502.2 is a general purpose method for the identification and simultaneous measurement of purgeable volatile organic compounds in finished drinking water, or drinking water in any treatment stage.

In the following application the use of Dani Master DHS/P&T operating in "purging mode" for the determination of VOCs **ACCORDING TO EPA METHOD 502.2** is presented. The Master DHS/P&T is the solution of choice to reach **HIGHER PRODUCTIVITY** and **INCREASED SENSITIVITY**.

*Do you wish to streamline the process, boost your productivity and save time and hassles?*

DANI Water Analyzer DWA-115 is the ready-to-go solution to attain the maximum performance in the shortest time for your analysis of Volatile Aromatic Organic Compounds in Drinking Water.

*Master* **DWA-115**

#### HIGHER PRODUCTIVITY

Overlapped thermostatted sample capability.  
Shorter baking phase.  
Use of disposable vials.

#### INCREASED SENSITIVITY

Superior sensitivity is obtained through the constant sweeping of the thermostatted sample, promoting the enrichment of the volatile compounds in the sorbent trap. The Master DHS/P&T offers a 100-fold increase in sensitivity over conventional headspace techniques.

#### THE CONFIGURATION MATCHES THE EPA METHOD 502.2

The system exceeds the low-level threshold required by the latest regulations.



### Analysis Conditions


#### Master GC Parameters:

Oven	35°C (8 min), 4°C/min, 240°C (1 min)
Detector	FID 250°C
Injector	SL/IN (220°C)
Carrier	Helium, 3,5 mL/min, (split 1:2)
Column	Vocol 60 m x 0.32 mm i.d. x 3µm d <sub>f</sub>

#### Master DHS/P&T Parameters - Purging Mode

Incubation	60°C
Stripping	3 min, 120 mL/min, Trap -10°C
Injection	3 min, Dew Stop 0°C, Trap 295°C
Baking	10 min, 80 mL/min Trap 300°C, Dew Stop 200°C
Transfer Line	250°C
Switching Valve	250°C
Trap Material	Tenax/Carbotrap/Carbosieve
Sample Volume	10mL

## Volatile Aromatic Organic Compounds in Drinking Water A Dynamic Approach to EPA Method 502.2



	COMPOUNDS	RSD%	MDL	R <sup>2</sup>	RANGE
			ppb		ppb
1	benzene	2.0	0.12	0.99952	0.4-10
2	toluene	2.1	0.12	0.99900	0.4-10
3	ethylbenzene	2.6	0.13	0.99851	0.4-10
4	<i>m</i> -xylene	7.9	0.09	0.99987	0.4-10
4	<i>p</i> -xylene	7.9	0.09	0.99987	0.4-10
5	<i>o</i> -xylene	4.0	0.15	0.99978	0.4-10
6	styrene	1.5	0.16	0.99935	0.4-10
7	<i>iso</i> -propylbenzene	1.7	0.25	0.99793	0.4-10
8	<i>n</i> -propylbenzene	2.4	0.25	0.99804	0.4-10
9	bromobenzene	0.9	0.32	0.99940	0.4-10
10	1,3,5-trimethylbenzene	0.3	0.19	0.99890	0.4-10
11	2-chlorotoluene	4.7	0.20	0.99932	0.4-10
12	4-chlorotoluene	1.0	0.20	0.99887	0.4-10
13	<i>tert</i> -butylbenzene	1.2	0.24	0.99894	0.4-10
14	1,2,4-trimethylbenzylene	0.5	0.18	0.99870	0.4-10
15	<i>sec</i> -butylbenzene	0.1	0.30	0.99829	0.4-10
16	<i>p</i> -isopropyltoluene	5.0	0.30	0.99815	0.4-10
17	1,3-dichlorobenzene	1.8	0.40	0.99879	0.4-10
18	1,4-dichlorobenzene	1.4	0.40	0.99889	0.4-10
19	<i>n</i> -butylbenzene	2.8	0.60	0.99674	0.4-10
20	1,2-dichlorobenzene	2.9	0.40	0.99874	0.4-10



# WATER ANALYSIS

## Dynamic Determination of Volatile Aromatic Hydrocarbons in Water

Application Note AN 119

Due to their toxicity and persistence in the environment, Volatile Aromatic Organic Compounds are particularly dangerous pollutants. About their possible effects on human health, it is important to highlight that some of them are mutagens, teratogens or carcinogens. For these reasons, government agencies require these contaminants to be monitored at progressive lower levels. The qualitative confirmation, quantitative accuracy and precision required in current regulations demand for high performing analytical solutions. In this work Master DHS/P&T coupled to Master GC demonstrates that the Purge&Trap is the technique of choice when in need to reach the minimum detectable levels required by law in force. This is demonstrated in the analysis reported below that shows **EXCELLENT CHROMATOGRAPHIC RESOLUTION** and **REPEATABILITY** with **NO RISK OF CROSS-CONTAMINATION**.

*Do you wish to streamline the process, boost your productivity and save time and hassles?*

DANI Water Analyzer DWA-119 is the ready-to-go solution to attain the maximum performance in the shortest time for your analysis of Volatile Aromatic Hydrocarbons in Water.

**Master DWA-119**

### EXCELLENT CHROMATOGRAPHIC RESOLUTION

The high desorption efficiency and the minimized sample path dead volume of the Master DHS/P&T guarantee unequalled chromatographic resolution.

### EXCEPTIONAL REPEATABILITY

A sophisticated control of all parameters and the highly precise electronic regulation of the purging gas flow-rate feature an unmatched repeatability (<2.5% RSD) and accuracy.

### NO RISK OF CROSS-CONTAMINATION

The Purge&Trap technique is based on the injection of vapors only, therefore ensuring a totally clean procedure and highly reliable results. The Master DHS/P&T, moreover, provides an automated cleaning cycle during the baking phase.



## Analysis Conditions

### Master GC Parameters:

Oven	40°C (3 min), 8°C/min, 160°C, 20°C/min, 230°C (5min)
Detector	FID 250°C
Injector	SL/IN (220°C)
Carrier	1.6 mL/min., (split 1:15)
Column	DN WAX 30m x 0.25 mm i.d. x 0.15 µm d <sub>f</sub>



### Master DHS/P&T Parameters - Purging Mode

Incubation	40°C
Stripping	15 min, 43 mL/min, Trap 30°C
Injection	1 min, Dew Stop 0°C, Trap 230°C
Baking	10 min, 80 mL/min Trap 280°C, Dew Stop 200°C
Transfer Line	200°C
Switching Valve	200°C
Trap Material	Tenax GR
Sample Volume	5 mL

## Dynamic Determination of Volatile Aromatic Hydrocarbons in Water

### Peak identification

- 1 benzene
- 2 toluene
- 3 ethylbenzene
- 4 *p*-xylene
- 5 *m*-xylene
- 6 *o*-xylene



Linearity was evaluated in the range 0.5-100 ppb

	benzene	toluene	ethylbenzene	<i>p</i> -xylene	<i>m</i> -xylene	<i>o</i> -xylene
MDL (ppb)	0.14	0.086	0.11	0.099	0.091	0.089

Assuming a Minimum Detectable Level 3 times the noise

	PEAK AREA					
	benzene	toluene	ethylbenzene	<i>p</i> -xylene	<i>m</i> -xylene	<i>o</i> -xylene
	53.87	63.38	49.48	56.12	61.81	62.43
	51.80	61.79	48.16	54.28	61.01	61.18
	53.30	62.72	48.82	54.70	60.91	61.94
	50.91	60.76	46.51	52.66	58.47	60.85
	51.43	59.54	46.73	53.07	59.01	60.17
	53.35	62.77	48.37	55.33	61.34	62.90
	53.08	63.07	49.42	55.97	61.26	62.69
	54.35	64.54	49.76	56.52	62.64	63.88
	52.13	62.57	48.38	55.08	61.40	61.75
Average	52.69	62.35	48.40	54.86	60.87	61.98
SD	1.17	1.48	1.15	1.33	1.32	1.14
RSD%	2.22	2.37	2.38	2.43	2.17	1.84

The repeatability obtained for a 50 ppb solution

# WATER ANALYSIS

## Automatic Dynamic Headspace and Purge&Trap Sampler for the Determination of VOCs in Water

### Application Note AN 120

Volatile Organic Compounds (VOCs) are organic chemicals for the most part dangerous for human health or harmful for the environment. Some of them are mutagens, teratogens or carcinogens. This is why VOCs are highly regulated by norms. The qualitative confirmation, quantitative accuracy and precision required in current regulations demand for high performing analytical solutions.

Purge&Trap is **THE PREFERRED TECHNIQUE** for its higher sensitivity, if compared to Static Headspace, when analyzing very low concentrations.

The aim of the following application is to show the **PRECISE RESULTS** and **THE HIGH SENSITIVE LEVELS OF ANALYSIS** obtained with the **EASY-TO-USE** Master DHS/P&T.

*Do you wish to streamline the process, boost your productivity and save time and hassles?*

DANI Water Analyzer DWA-120 is the ready-to-go solution to attain the maximum performance in the shortest time for your analysis of VOCs in Water.

**Master DWA-120**

#### THE PREFERRED TECHNIQUE FOR THE ANALYSIS OF VOCs

The Master DHS/P&T provides the highest sensitivity and the capability to analyze water samples at very low concentrations using the Purge&Trap technique.

#### ACCURATE, PRECISE, AND EASY-TO-USE SOLUTION

Unlike conventional Purge&Trap systems, the overlapped sample thermostating capability, the shorter baking phase, and the use of disposable vials allow the Master DHS/P&T to maximize productivity with simplicity.

#### HIGH SENSITIVITY

Highest sensitivity and excellent performances are guaranteed by the use of the ECD - Electron Capture Detector.



### Analysis Conditions

#### Master GC Parameters:


Oven	40°C (6 min), 14°C/min., 230°C (11 min)
Detector	FID, 250°C, ECD 300°C
Carrier	Helium, 1 mL/min, (split 1:10)
Column	Rtx-VMS, 60 m x 0.25mm i.d. x 1.4 µm d <sub>f</sub>
Injector	SL/IN (230°C)

#### Master DHS/P&T Parameters - Purging Mode

Incubation	40°C
Stripping	10 min, 60 mL/min, Trap 20°C
Injection	1 min, Dew Stop 0°C, Trap 230°C
Baking	10 min, 80 mL/min Trap 250°C, Dew Stop 200°C
Transfer Line	200°C
Switching Valve	200°C
Trap Material	Tenax/Carbotrap/Carbosieve
Sample Volume	10 mL




## Automatic Dynamic Headspace and Purge&Trap Sampler for the Determination of VOCs in Water



The 50 ppb standard solution was analyzed with FID

Peak identification	
1	1,1-dichlorethene
2	methylene chloride
3	<i>trans</i> -1,2-dichloroethene
4	<i>cis</i> -1,2-dichloroethene
5	chloroform
6	carbon tetrachloride
7	1,1,1-trichloroethene
8	benzene
9	1,2-dichloroethane
10	trichloroethene
11	1,2-dichloropropene
12	bromodichloromethane
13	<i>cis</i> -1,3-dichloropropene
14	toluene
15	<i>trans</i> -1,3-dichloropropene
16	tetrachloroethane
17	1,1,2-trichloroethane
18	dibromochloromethane
19	<i>m</i> -xylene
20	<i>p</i> -xylene
21	<i>o</i> -xylene
22	bromoform
23	1,4-dichlorobenzene



The 0.01 ppb standard solution was analyzed with ECD

# WATER ANALYSIS

## Determination of BTEX and Styrene in Water Using Static Headspace

### Application Note AN 049

BTEX and Styrene are Volatile Organic Compounds derived from petroleum and from the emissions of motor vehicles. These compounds are known for the contamination of soils and groundwater and for their harmful effects on human health. Because of their volatility, Static Headspace is an advisable technique for the analysis of BTEX and Styrene in water. Moreover, unlike other sample handling methodologies, Static Headspace is an easy-to-use, solvent-free, and robust technique. The aim of the following work is to show that by using the state-of-the-art DANI Master SHS, it is possible to achieve **HIGHLY ACCURATE RESULTS** at ppb levels even for small sample volumes at low concentrations. This approach **ELIMINATES THE RISKS OF CARRY-OVER AND CROSS-CONTAMINATION** and allows to save time **INCREASING LABORATORY PRODUCTIVITY**.

*Do you wish to streamline the process, boost your productivity and save time and hassles?*

DANI Water Analyzer DWA-049 is the ready-to-go solution to attain the maximum performance in the shortest time for your analysis of BTEX and Styrene in Water.

**Master DWA-049**

#### INCREASED LABORATORY PRODUCTIVITY

120 Sample Tray guarantees the highest sample capacity for the highest productivity.

#### SAMPLE INTEGRITY PRESERVATION

The sample flow path of the Master SHS is entirely chemically inert and can be thermostatted to high temperatures. These features eliminate analytical carryover and maintain sample integrity.

#### HIGH ACCURACY AT VERY LOW CONCENTRATION

The Valve&Loop technique is the most reliable and used technique which is capable of highly repeatable results. Master SHS can guarantee outstanding repeatability and avoid the risk of false results, sample loss or recondensation.



### Analysis Conditions


#### Master GC Parameters:

Oven	40°C, 5°C/min., 110°C (5 min), 20°C/min., 200°C
Detector	FID, 250°C
Injector	SL/IN, 150°C
Injection mode	split, split flow 10 mL/min, split ratio 1:1
Carrier	Helium, 10 mL/min
Column	DN-WAX, 25m x 0.53 mm i.d. x 1.2 µm d <sub>f</sub>
Sample Volume	BTEX and Styrene in water 10 mL

#### MasterSHS parameters:

Manifold	85°C
Oven	75°C
Transfer Line	85°C
Incubation Time	30 min
Aux. Gas	0.7 bar


## Determination of BTEX and Styrene in Water Using Static Headspace



Chromatogram of a 10 µg/L standard solution

	benzene	toluene	ethylbenzene	p-xylene	m-xylene	o-xylene	styrene
1	14.21	16.63	17.59	17.56	18.05	16.82	14.99
2	14.27	16.40	17.44	17.33	17.74	16.79	14.87
3	14.21	16.16	17.12	16.94	17.34	16.54	14.63
4	14.12	15.87	17.05	16.82	17.17	16.29	14.58
5	14.01	16.04	17.24	16.84	17.34	16.36	14.85
6	13.87	16.05	17.02	16.87	17.17	16.32	14.71
Average	14.12	16.19	17.24	17.06	17.48	16.52	14.77
SD	0.15	0.28	0.23	0.31	0.34	0.24	0.16
RSD%	1.06	1.72	1.32	1.81	1.97	1.43	1.05

Repeatability and RSD%, obtained for a 10 µg/L standard solution, are calculated on six repetitions



System linearity calculated in a range from 0.2 to 100 µg/L



# WATER ANALYSIS

## Phenols: FAST GC ANALYSIS

### Application Note AN 007

Phenols are a class of very common chemical compounds. They can be found in the natural world and they are also used as raw materials and additives for industrial purposes in preservatives, insecticides, and plastics. Releases of phenols in water result from wastewater from manufacturing industries and from commercial use of phenol and phenol-containing products. Phenols have been detected in surface waters, groundwater, drinking water and at hazardous waste sites.

They represent a danger to the environment and to human health. In fact, phenols are hematotoxic and hepatotoxic, provoke mutagenesis and carcinogenesis towards humans and other living organisms. The presence of phenols should be limited to 0.3 milligrams per liter of water to protect human health from the possible harmful effects of exposure to phenol by drinking water and/or eating contaminated water plants and animals.

For this reason U.S. EPA takes into account the analysis of phenols in a variety of methods including EPA Method 604 and 8041.

The following analysis demonstrates a **FAST GC ANALYSIS** for eleven target compounds in less than five minutes showing a **COST-EFFECTIVE METHOD** with **HIGH RESOLUTION POWER** and **EXCELLENT ACCURACY**.

#### COST-EFFECTIVE METHOD FOR THE FAST DETERMINATION OF PHENOLS IN WATER

Conventional GC average analysis time: 30 minutes.

DANI Master GC analysis time : less than 5 minutes.

#### HIGH RESOLUTION POWER AND EXCELLENT ACCURACY

The Fast Dedicated Column with narrower internal diameter and thinner stationary phase films features faster analysis time while maintaining proper resolving power.

#### SUPERIOR RETENTION TIME STANDARD DEVIATION

Great precision is obtained with an average Retention Time Standard Deviation of 0.0013 min




### Analysis Conditions

#### Master GC Parameters:

Oven	80°C (1min) 40°C/min, 220°C, (0.5min)
Detector	FID 400° C
Injector	PTV 80°C, 600°C/min, 400°C
Split Flow	25mL/min, split ratio 1:50
Carrier	H <sub>2</sub> 0.5mL/min
Column	DN 5 FAST 15m x 0.10mm i.d. x 0.10 µm d <sub>i</sub>
Sample Volume	0.5 µL

## Phenols: FAST GC ANALYSIS



Compounds	Retention Time Std dev. (min)	Compounds	Retention Time Std dev. (min)
phenol	0,0013	2,4,6-trichlorophenol	0,0012
2-chlorophenol	0,0013	2,4-dinitrophenol	0,0016
2-nitrophenol	0,0016	4-nitrophenol	0,0010
2,4-dimethylphenol	0,0015	2-methyl-4-dinitrophenol	0,0013
2,4-dichlorophenol	0,0013	pentachlorophenol	0,0015
4-chloro-3-methylphenol	0,0012		

# WATER ANALYSIS

## Poly Aromatic Hydrocarbons (PAHs) - Fast GC Application

### Application Note AN 003

PAHs are by-products of petroleum processing or combustion. Many of these compounds are highly carcinogenic and organic pollutants at relatively low levels. Although they are nearly insoluble in water, their highly hazardous nature justifies the need for monitoring their presence in potable waters and wastewaters. As proof of this, the Environmental Protection Agency has included 16 PAHs on its list of priority pollutants to be monitored.

U.S. EPA 8100 method provides gas chromatographic conditions for the detection of ppb levels of certain polyaromatic hydrocarbons.

The aim of the following application is to present the **FAST ANALYSIS** of 18 representatives of the PAHs class of compounds at concentrations that **MEET THE EPA METHOD REQUIREMENTS** in less than 10 minutes. The results show **OUTSTANDING RESOLUTION POWER**.

#### PAHs FASTEST ANALYSIS

Fast Dedicated Column (DN-PAH-FAST).  
Fast Acquisition Rate of the Detector (300 Hz).

#### OUTSTANDING RESOLUTION POWER

The system, along with the optimal control of the oven temperature and the fast detector, assures an outstanding resolution power, unprecedented for all chromatographic measurements.

#### REGULATORY COMPLIANT RESULTS




### Analysis Conditions

#### Master GC Parameters:


Oven	140°C (0.5min) 30°C/min, 220°C, 15°C/min, 300°
Detector	FID 400° C
Injector	PTV 80°C, 600°C/min, 400°C
Split Flow	50 mL/min, split ratio 1:100
Carrier	H <sub>2</sub> 0.5mL/min
Column	DN PAH FAST 15m x 0.10mm i.d. x 0.10 µm d <sub>f</sub>
Sample Volume	0.5 µL



## Poly Aromatic Hydrocarbons (PAHs) - Fast GC Application



phenanthrene and anthracene resolution value : 1.35



# WATER ANALYSIS

## Fast GC approach for PCBs Determination

### Application Note AN 121

Polychlorinated biphenyls are a class of organic compounds known for their high level of toxicity and classified as persistent organic pollutants. Thanks to their useful characteristics such as non-flammability, heat resistance, insulation and chemical stability, in the past they were extensively used as coolants and dielectric fluids, stabilizing additives in PVC and plastic products, reactive flame retardants, sealants, paints, etc.

The toxicity associated to PCBs was recognized and known very soon, since before their first commercial production in 1970s. Nevertheless, PCB production was banned by the United States Congress only in 1979 and by the Stockholm Convention on Persistent Organic Pollutants in 2001.

The **EPA METHOD 8082** is used to determine the concentration of PCBs in extracts from solid and aqueous matrices. The analytical protocol is based on conventional gas chromatography coupled to electron capture detection technique. This method generally requires 20-30 minutes for the chromatographic separation of these compounds. Fast gas chromatography, typically involving 100  $\mu\text{m}$  i.d. and 10 m columns, represents a powerful alternative to conventional GC, allowing to achieve equivalent **RESOLUTION IN SIGNIFICANTLY SHORTER ANALYSIS TIME**.

In this application, fast GC is applied to the analysis of PCBs. The technique is implemented on DANI Master GC, which operates with short narrow bore columns and fast ECD to guarantee high resolution and sensitivity. The data obtained confirm the suitability of the technique for the routine analysis of this kind of compounds.

*Do you wish to streamline the process, boost your productivity and save time and hassles?*

DANI Water Analyzer DWA-121 is the ready-to-go solution to attain the maximum performance in the shortest time for your analysis of PCBs in Water.


**Master DWA-121**

#### **SIGNIFICANT REDUCTION OF THE ANALYSIS TIME WITHOUT LOSS OF RESOLUTION**

High acquisition rate up to 300 Hz  
High separation power of the column

#### **FAST AND CONVENTIONAL ANALYSIS IN A UNIQUE SYSTEM**

The versatile and flexible Master GC is uniquely designed to perform both conventional and fast gaschromatographic analyses.



#### **ACHIEVEMENT OF THE LIMITS SET BY THE EPA METHOD**

EPA Method 8082 mandates quantitation limits down to 0,17 ng/L

	Conventional GC		Fast GC
Column	DN 5-30m x 0.25mm i.d. x 0.25 $\mu\text{m}$ d <sub>f</sub>	DN 5 - 10m x 0.1mm i.d. x 0.2 $\mu\text{m}$ d <sub>f</sub>	DN 5 - 5m x 0.1mm i.d. x 0.2 $\mu\text{m}$ d <sub>f</sub>
PTV Injector	50°C, 600°C/min, 320°C (2min)	50°C, 600°C/min, 320°C (2min)	50°C, 600°C/min., 320°C (2min.)
Oven	120°C, 10°C/min, 300°C (4min)	120°C, 25°C/min, 200°C, 20°C/min., 300°C/min (2 min)	120°C, 25°C/min., 200°C, 20°C/min., 300°C/min. (2 min.)
Carrier Gas (Helium) Flow Rate	1mL/min	0.5 mL/min	0.5 mL/min.
Split Ratio	1:10	1:50	1:50
ECD Detector	320°C	320°C	320°C
Digital Aquisition Rate	25 Hz	300 Hz	300 Hz
Injection Volume	1 $\mu\text{L}$	0.5 $\mu\text{L}$	0.5 $\mu\text{L}$


## Fast GC approach for PCBs Determination

Compound	Repeatability (min)		Water Sample Limit ng/L	
	RT (SD)	Area (RSD)	LOD	LOQ
1) 2,4,5,6-tetrachloro- <i>m</i> -xylene (S.S)	0.001	1.011	0.052	0.175
2) 2,3-dichlorobiphenyl	0.001	0.819	3.030	10.101
3) 2,2',5'-trichlorobiphenyl	0.000	1.021	5.607	18.692
4) 2,4',5-trichlorobiphenyl	0.002	1.267	4.511	15.038
5) 2,2',5,5'-tetrachlorobiphenyl	0.001	1.667	3.593	11.976
6) 2,2',3,5'-tetrachlorobiphenyl	0.001	1.584	2.439	8.130
7) 2,3',4,4'-tetrachlorobiphenyl	0.002	1.101	2.120	7.067
8) 2,2',4,5,5'-pentachlorobiphenyl	0.001	1.757	2.120	7.067
9) 2,2',3,4,5'-pentachlorobiphenyl	0.002	1.319	1.354	4.515
10) 2,3,3',4',6-pentachlorobiphenyl	0.002	0.942	1.304	4.348
11) 2,2',3,5,5',6-hexachlorobiphenyl	0.002	1.509	1.017	3.390
12) 2,2',4,4',5,5'-hexachlorobiphenyl	0.002	1.770	0.789	2.632
13) 2,2',3,4,5,5'-hexachlorobiphenyl	0.001	1.572	0.589	1.963
14) 2,2',3,4,4',5'-hexachlorobiphenyl	0.001	1.274	0.913	3.044
15) 2,2',3,4,5,5',6-heptachlorobiphenyl	0.002	1.725	0.557	1.857
16) 2,2',3,4,4',5',6-heptachlorobiphenyl	0.002	0.917	0.507	1.691
17) 2,2',3,4,4',5,5'-heptachlorobiphenyl	0.002	0.950	0.489	1.630
18) 2,2',3,3',4,4',5-heptachlorobiphenyl	0.001	1.385	0.503	1.675
19) 2,2',3,3',4,4',5,5'-nonachlorobiphenyl	0.002	1.132	0.500	1.668
20) decachlorobiphenyl (I.S.)	0.002	1.364	0.039	0.130

Repeatability, LOD and LOQ calculated for each target compound were obtained with the 10 m column (comparable results were achieved with the 5 m column)

**CONVENTIONAL GC**  
30 m x 0.25 mm column  
acq. rate 25 Hz

$P_{W05}=0.040$  min




$S/N=161.128$

**FAST GC**


10 m x 0.1 mm column  
acq. rate 100 Hz

$P_{W05}=0.023$  min



$S/N=480.275$

Peak width and S/N comparison : Fast GC provides narrower peaks and a higher S/N ratio



Analysis of a 20 PCB congeners mixture (500 ppb): analysis time and resolutions obtained with three different columns



# WATER ANALYSIS

## Fuel Oxygenates in Water

### Application Note AN 122

Fuel oxygenates, primarily ethers and alcohols, are added to gasoline to enhance the octane content and to improve air quality reducing the emission of pollutants, particularly carbon monoxide.

Fuel oxygenates can be found in aquifers as contaminants. Their introduction to the environment can be accidental through the release from underground pipelines, tanks and gasoline spills. Industrial wastewater as well as the petroleum fuel cycle process can also represent sources of oxygenates to the water.

The proof of the presence of oxygenates in drinking water has raised serious concern regarding the taste and odor aspects. Concerns also raised about possible human-health implications.

In the following application **INCREASED ANALYTICAL PERFORMANCES** and **HIGH RELIABILITY** are obtained through the use of the Master GC coupled to the Master DHS/P&T.

*Do you wish to streamline the process, boost your productivity and save time and hassles?*

DANI Water Analyzer DWA-122 is the ready-to-go solution to attain the maximum performance in the shortest time for your analysis of Fuel Oxygenates in Water.

**Master DWA-122**

#### INCREASED ANALYTICAL PERFORMANCES

Master DHS/P&T incorporates the innovative and ingenious Dew Stop device which efficiently removes water regardless of the analytes, maintaining volatile compounds recovery unaffected.

#### HIGH RELIABILITY

The minimal sample handling required by the system, along with the complete automation of all process steps, ensures highly reliable and reproducible results minimizing operators errors.



## Analysis Conditions

### Master GC Parameters:


Oven	35°C , 6°C/min., 90°C , 45°C/min., 210°C (15 min.)
Detector	FID 250°C
Injector	SL/IN 200°C
Carrier	Helium, 1.2 mL/min., (split 1:20)
Column	Vocol 60 m x 0.25 mm i.d. x 1.5 µm d <sub>f</sub>

### Master DHS/P&T Parameters - Purging Mode


Incubation	0.5 min., 60°C, Shaking Fast
Stripping	11 min., 40 mL/min., Trap 0°C
Injection	2 min., Dew Stop 0°C, Trap 300°C
Baking	5 min., 150 mL/min. Trap 310°C, Dew Stop 200°C
Transfer Line	170°C
Switching Valve	170°C
Trap Material	Carbopack B/Carboxen 1000
Sample Volume	10mL



## Fuel Oxygenates in Water



Chromatogram of a standard solution



System linearity calculated in a range from 0.02 to 200 ppb

# WATER ANALYSIS

## Hydrocarbon Oil Index in Water

### ISO 9377-2 method

#### Application Note AN 123

The determination of the Hydrocarbon Oil Index is mandatory for the environment and human health protection. The ISO 9377-2 is the official European method for oil and grease determination in water. This test is a gas chromatographic method suitable for surface water, wastewater and water from sewage treatment. The goal of the following application is to show a system configuration to fulfill the requirements of the method **EASILY** and **IN A VERY SHORT TIME**. Fast GC analysis is demonstrated to be **A RELIABLE, PROVEN AND AUTOMATED TECHNIQUE** able to **IMPROVE LABORATORY PRODUCTIVITY**. This is why the proposed system configuration is the perfect solution for those laboratories that are constantly faced with the need to maximize sample throughput without sacrificing the accuracy of the results.

*Do you wish to streamline the process, boost your productivity and save time and hassles?*

DANI Water Analyzer DWA-123 is the ready-to-go solution to attain the maximum performance in the shortest time for your analysis of the Hydrocarbon Oil Index in Water.

**Master DWA-123**

#### RELIABLE, PROVEN AND AUTOMATED TECHNIQUE

The Master GC, uniquely designed to perform both conventional and fast gas chromatographic analyses, deliver unsurpassed analytical capabilities. Moreover, unlike other commercially available fast gas chromatographs, the Master GC offers guided diagnostic and maintenance procedures supporting the user in the preservation of the system precision.

#### IMPROVED LABORATORY PRODUCTIVITY

Fast analysis time, fast results and the consequent reduction of the cost per analysis are the driving factors for every environmental laboratory.

#### COMPLIANT WITH ISO 9377-2 METHOD

The proposed configuration fulfills all the requirements of the method.




## Analysis Conditions

### Master GC Parameters:


Oven	125°C, 70°C/min., 175°C, 50°C/min., 300°C, 35°C/min., 350°C (1 min.)
Detector	FID 380°C
Injector	PTV 100°C, 999°C/min., 380°C (2 min.)
Split Flow	1:30
Carrier	He 1,89 bar
Column	DN-5 FAST 5 m x 0,1 mm i.d. x 0,1 µm d <sub>f</sub>
Sample Volume	1 µL

## Hydrocarbon Oil Index in Water ISO 9377-2 method


### Alkanes C7 - C40



### Alkanes C8 - C40



### Alkanes C10 - C40



# WATER ANALYSIS

## Determination of Phtalates

### Application Note AN 052

Phtalate esters are contaminants mainly used in a large variety of products such as children toys, entering coatings of pharmaceutical pills, cosmetics, detergents, film formers and, more generally, plastic products. Recent studies link phtalates to different human diseases, from disruption for the endocrine system to cancer.

Phtalates are easily released into the environment due to the plastic breakdown and aging. Due to their massive presence in the environment, phtalates are also commonly found in groundwater.

Direct or indirect exposure to these compounds may cause health diseases. Phtalates can be found almost everywhere; for this reason EPA has developed the method 606 in order to quantify them. EPA Method 606 is a gas chromatographic method applicable to the determination of phtalate esters in municipal and industrial discharges.

The application below shows a **SIMPLE** and **RELIABLE** solution for the analysis of Phtalates.

#### **SIMPLE AND RELIABLE DETERMINATION OF PHTALATES**

Unparalleled and reliable chromatographic accuracy and precision are guaranteed by the patented Digital Flow Control. All the parameters can be easily set up and controlled by an intuitive touchscreen.

#### **ONLY 10 MINUTES ANALYSIS TIME**

The short analysis run times and higher performances of the Master GC significantly reduce laboratory operating costs.




## Analysis Conditions

### Master GC Parameters:

Oven	150°C - 15°C/min - 270°C
Detector	FID 300°C
Injector	SL/IN 250°C
Injection Mode	Splitless
Carrier	He 20 ml/min
Column	DN-1 15m x 0,53mm i.d. x 1,50µm d <sub>f</sub>
Sample Volume	1.0 µL



Determination of Phthalates  
EPA Method 606



Peak Identification	
1	Dimethyl-phthalate
2	Diethyl-phthalate
3	Di-n-butyl phthalate
4	Butyl benzyl phthalate
5	Bis (2-ethylhexyl) phthalate
6	Di-n-octyl phthalate

# WATER ANALYSIS

## Determination of Haloethers

### EPA Method 611/8111

#### Application Note AN 069

In the perspective of a study about all identifiable effects on health and welfare which may be expected from the presence of pollutants in any body of water, including ground water, a special attention is reserved to haloethers.

Haloethers are pollutant compounds mostly manufactured and they are used as solvents, chemical intermediates, soil fumigants, pesticides, fungicides, etc. Moreover, haloethers are characterized by their persistence in natural surface waters and can be adsorbed by organic-rich sediments and bioaccumulated in fish.

For the above-mentioned reasons haloethers are under investigation as a possible cause for different types of human diseases. EPA method 8111 provides gas chromatographic conditions for the detection of ppb concentration of haloethers in water and soil or ppm concentration in waste samples.

**METHOD REQUIREMENTS ARE ACHIEVED** in the following analysis with **MAXIMUM PRECISION AND ACCURACY**.

#### MAXIMUM PRECISION

The patented DFC - Digital Flow Control automatically adjusts the carrier gas flow to compensate the ambient temperature and pressure providing constant retention time, enhanced repeatability and extreme precision.

#### EXTREMELY FLEXIBLE SYSTEM

The advanced modular design of the Master GC components features outstanding flexibility and upgradeability. Any GC configuration can be easily modified and/or upgraded.

#### UNPARALLELED CHROMATOGRAPHIC ACCURACY

The patented DFC and the optimal control of the oven temperature assure outstanding retention time repeatability, unprecedented for all chromatographic measurements.




## Analysis Conditions

### Master GC Parameters:

Master GC Oven	100°C - 15°C/min - 300°C
Detector	FID 300°C
Injector	SL/IN 250°C
Split Ratio	1:50
Carrier	He 2,90 psi
Column	DN-5 15m x 0,53mm i.d x 1,50µm d <sub>f</sub>
Sample Volume	0.2 µL

Determination of Haloethers  
EPA Method 611/8111



**Peak identification**

- |   |                                      |
|---|--------------------------------------|
| 1 | <i>bis</i> (2-chloroethyl) ether     |
| 2 | <i>bis</i> (2-chloroisopropyl) ether |
| 3 | <i>bis</i> (2-chloroethoxy)methane   |
| 4 | 4-chlorophenylphenyl ether           |
| 5 | 4-bromophenyl phenyl ether           |



# WATER ANALYSIS

## Determination of Nitrosamines

### EPA Method 607

#### Application Note AN 067

Nitrosamines are a family of compounds used in the manufacture of rubber, cosmetics, pesticides, leather, etc. They can be also found in tobacco, cured meats, and beer.

Thanks to the massive studies involving Nitrosamines, it has been demonstrated that these compounds are mutagens and carcinogens. In regard to this, US Environmental Protection Agency and worldwide environmental and health related government agencies have imposed restrictions on the use of these substances. It is therefore essential to be able to count on a reliable and accurate solution that ensures **EXCELLENT RESULTS IN TERMS OF REPRODUCIBILITY AND PEAK SEPARATION**.

EPA Method 607 is a gas chromatographic method applicable to the determination of certain nitrosamines in municipal and industrial discharges. The following work shows excellent results that **MATCH THE REQUIREMENTS OF THE METHOD**.

#### EXCELLENT RESULTS IN TERMS OF REPRODUCIBILITY AND PEAK SEPARATION

The column combined with the fast GC oven temperature generates sharper peaks. The proprietary PTV Injector achieves extremely fast heating rates and rapid cool down with ambient air. After injection, the PTV can be programmed to decrease the split flow and save carrier gas. The DFC, in addition, adjusts the carrier gas flow providing ambient temperature and pressure compensation. These features provide constant retention time and unmatched reproducibility.



#### RESULTS MEETING EPA METHOD REQUIREMENTS

The proposed configuration provides reliable and repeatable results in compliance with the EPA method 607.


#### Analysis Conditions

##### Master GC Parameters:

Oven	40°C - 20°C/min - 240°C
Detector	FID 280°C
Injector	PTV 50°C, 600°C/min, 240°C
Carrier	H2 10 ml/min
Column	DN-5 15m x 0,53mm i.d x 1,50 µm d <sub>i</sub>
Sample Volume	1.0 µL



Determination of Nitrosamines  
EPA Method 607

**Peak identification**

- |   |  |
|---|--|
| 1 | <i>n</i> -nitrosodimethylamine             |
| 2 | <i>n</i> -nitrosodi- <i>n</i> -propylamine |
| 3 | <i>n</i> -nitrosodiphenylamine             |

# WATER ANALYSIS

## Determination of Phenols and Chlorophenols

### EPA Method 604

#### Application Note AN 066

Phenols exist in the environment as products of the chemical, petrol, tinctural and pharmaceutical industries and as a consequence of a number of pesticides and the generation of industrial sewages.

Phenols are one of the first compounds reported into the List of Priority Pollutants by the US Environmental Protection Agency for their toxicity.

EPA Method 604 is a flame ionization detector gas chromatographic (FIDGC) method for the determination of phenols and certain substituted phenols in municipal and industrial discharges.

The analysis here below shows the easy **ACHIEVEMENT OF THE EPA METHOD REQUIREMENTS**, thanks to a **FLEXIBLE SYSTEM, SIMPLE TO SET UP AND TO CONTROL**.

#### **FLEXIBILITY**

The advanced modular design of the Master GC components features outstanding flexibility and upgradeability. Any GC configuration can be easily modified or upgraded. The Master GC allows the assembly of up to three injector and three detectors simultaneously.

#### **QUICK AND EASY SET UP**

The Master GC incorporates an intuitive and easy-to-use touchscreen interface that provides quick and easy set up and control. The system can also be controlled by the functional and user-friendly CLARITY™ Chromatography Station.

#### **COMPLIANCE WITH REGULATORY NORMS**

The proposed configuration provides reliable and repeatable results in compliance with the EPA method 604.




## Analysis Conditions

### Master GC Parameters:

Oven	110°C - 8°C/min - 280°C
Detector	FID 300°C
Injector	SL/IN 300°C
Carrier	H <sub>2</sub> 8,70 psi
Split Flow	1:100
Column	DN-5 25m x 0.32mm i.d. x 1,00 µm d <sub>f</sub>
Volume Injected	1.0 µL

Determination of Phenols and Chlorophenols  
EPA Method 604



**Peak identification**

1	phenol	7	2,4,6-trichlorophenol
2	2-chlorophenol	8	2,4-dinitrophenol
3	2-nitrophenol	9	4-nitrophenol
4	2,4-dimethylphenol	10	2-methyl-4,6-ditrophenol
5	2,4-dichlorophenol	11	pentachlorophenol
6	4-chloro-3-methylphenol		



# WATER ANALYSIS

## Determination of Organochlorinated Pesticides

### EPA Method 608/8081

#### Application Note AN 063

Organochlorinated pesticides have a long history of widespread use and are persistent organic pollutants. Traces of these pesticides can still be found in the environment in the top layer soils after more than twenty years they have been banned. They have significant toxicity to plants, animals and humans, accumulating in food chains. It is therefore important to rely on an **ACCURATE CHROMATOGRAPHIC SEPARATION** and to obtain an **EXACT QUANTIFICATION** easily even in complex matrices.

EPA Method 608 is a gas chromatographic (GC) method applicable to the determination of certain organochlorinated pesticides and PCBs in municipal and industrial wastes.

EPA Method 8081 is used to determine the concentrations of various organochlorinated pesticides in extracts from solid and liquid matrices.

The **REQUIREMENTS OF BOTH METHODS ARE ACHIEVED** in the following analysis.

#### ACCURATE SEPARATION FOR PESTICIDES WITHOUT RISK OF DISCRIMINATION AND DEGRADATION

Unique characteristics of DANI PTV : the sample is introduced by cold injection followed by vaporization eliminating possible discrimination or degradation.

#### EASY QUANTIFICATION OF ORGANOCHLORINATED PESTICIDES

Selective Detector (ECD) offers excellent performances in the determination of pesticides. The sensitivity of the ECD enables it to provide unmatched performances for tough applications.



#### COMPLIANT WITH EPA METHODS 608/8081

The proposed configuration provides reliable and repeatable results in compliance with the EPA method 608/8081.


## Analysis Conditions

### Master GC Parameters:

Oven	65°C - 20°C/min - 150°C - 7°C/min - 260°C
Detector	ECD 280°C
Injector	PTV 50°C, 600°C/min, 260°C
Carrier	H <sub>2</sub> 8,70 psi
Column	DN-5 25m x 0,32mm i.d. x 0,25 µm d <sub>f</sub>
Sample Volume	1.0 µL



## Determination of Organichlorinated Pesticides EPA Method 608/8081



Peak identification			
1	$\alpha$ - BHC	9	4,4' DDE
2	$\beta$ - BHC	10	dieldrin
3	$\gamma$ - BHC	11	endrin
4	$\delta$ - BHC	12	4,4' DDD
5	heptachlor	13	endosulfan II
6	aldrin	14	endrin aldehyde
7	heptachlor epoxide	15	4,4' DDT
8	endosulfan I	16	endosulfan sulfate

# WATER ANALYSIS

## Determination of Organophosphorus Pesticides

### Application Note AN 095

Organophosphorus Pesticides are among the most widely used class of pesticides thanks to their high efficacy against pests. They are also well known for their poisoning effects on human health as the over-exposure to organophosphorus pesticides may cause irreversible damage to the nervous system and have neurotoxic effects on developing organisms. They can be, in fact, absorbed by inhalation, ingestion, and dermal absorption.

These compounds represent a concrete risk for the environment, also. Industrial waste, seepage from buried toxic wastes, and contamination during spraying operations, they all can be considered as possible ways for their introduction into the water.

EPA Method 622 is a gaschromatographic (GC) method applicable to the determination of certain organophosphorus pesticides in industrial and municipal discharges as provided under 40 CFR 136.1.

The analysis below is an example, applicable to complex mixtures, that shows how to **REACH THE METHOD DETECTION LIMITS**. These **OUTSTANDING RESULTS** can be achieved thanks to an **INCREASED PEAK RESOLUTION** and the use of a particularly selective and sensitive detector.

#### INCREASED PEAK RESOLUTION

The column (DN-68) is dedicated to the analysis of phosphorus pesticides and generates sharper peaks which result in higher signal and greater signal-to-noise ratios.

#### OUTSTANDING RESULTS IN TERMS OF S/N RATIO EVEN FOR NOT COMPLETELY PURIFIED SAMPLES

The Flame Photometric Detector selectively detects compounds containing sulfur or phosphorus. It is the detector of choice for the detection of phosphorus pesticides by virtue of its selectivity and sensitivity.

#### THE CONFIGURATION MATCHES THE EPA METHOD 622

The proposed configuration provides reliable and repeatable results in compliance with the EPA method 622.




## Analysis Conditions

### Master GC Parameters:

Master GC Oven	100°C (1 min) - 5.5°C/min - 230°C - 30°C/min - 270°C
Detector	FPD 140°C
Injector	SL/IN 270°C
Injection Mode	Splitless
Carrier	H <sub>2</sub> 10,15 psi
Column	DN-68 25m x 0.32mm i.d. x 0.25 µm d <sub>f</sub>
Sample Volume	1.0 µL

## Determination of Organophosphorus Pesticides



Peak identification	
1	Trichlorfon
2	Phosdrin
3	Tionazine
4	Ethoprophos
5	Phorate
6	Phonophos
7	Diazinone
8	CH <sub>3</sub> -Chlorpyrifos
9	CH <sub>3</sub> -Parathion
10	CH <sub>3</sub> -Pirimiphos
11	Chlorpyrifos
12	Malathion
13	Parathion
14	Pirimiphos
15	Quinalphos
16	C <sub>2</sub> H <sub>5</sub> -Bromophos
17	Metidathion
18	Ethion
19	Trithion
20	CH <sub>3</sub> -Azinphos
21	C <sub>2</sub> H <sub>5</sub> -Azinphos