

Chloramphenicol ELISA kit

Chloramphenicol is primarily bacteriostatic. It inhibits bacterial protein synthesis.

The drug is associated with random cases of aplastic anaemia, a serious blood disorder for which there is no cure, and which is usually fatal. In addition, there are concerns related to potential carcinogenicity and genotoxicity of the drug as well as the potential to cause antimicrobial resistance. therefore this drug is banned for use in food-producing animals in the majority of countries. According to a Codex Alimentarius decision, CAP should be banned from food production.

Assay principle

Competitive enzyme immunoassay for the quantitative analysis of chloramphenicol (CAP). This CAP-enzyme immunoassay is suitable for the detection of CAP and its glucuronide (major metabolite in urine).

Applications

Urine, milk, honey, serum, plasma, feed, eggs, tissue, shrimps, fish and water.

Standard curve range

0.1-2 ng/ml

Detection limits

Serum, water and urine: 0.1 ppb

Shrimps and fish: 0.02 ppb

Tissue, feed and egg: 0.05 ppb

Milk: 0.01 ppb

Honey: 0.025 ppb

Measuring range

Serum, water and urine: 0.1 – 2 ppb

Shrimps and fish: 0.02 – 0.4 ppb

Tissue, feed and egg: 0.05 – 1 ppb

Milk: 0.01 – 0.2 ppb

Honey: 0.025 – 0.5 ppb

Specificity (% cross-reactivity)

Chloramphenicol 100%

Chloramphenicol glucuronide 70%

Sample preparation

Serum and plasma: extraction with solvents.

Urine:

method I: adjustment the pH if necessary.

method II: enzymatic hydrolysis, extraction.

Tissue, feed, shrimps, fish and egg: homogenization, extraction with solvents and centrifugation

Honey: dilution, extraction with solvents and centrifugation.


Milk: centrifugation, extraction with solvent and dilution.

Water: adjustment the pH if necessary.

Assay time

45 minutes, sample preparation not included.

Example of standard curve



Kit contents

- Microtiter 96 well plate, 12 strips with 8 breakable wells)
- Chloramphenicol standard solutions: 0, 0.1, 0.2, 0.5, 1, 2 ng/ml, ready to use
- Enzyme conjugate, concentrated
- Enzyme conjugate diluent, ready to use
- Anti-chloramphenicol antibody, ready to use
- Dilution buffer, concentrated
- Developing solution, ready to use
- Washing buffer, concentrated
- Stop solution, ready to use
- CAP spiking solution
- Instruction

Materials not provided

- Distilled water
- Ethylacetate
- Acetic acid (urine: method II)
- β -glucuronidase from *Helix pomatia* (urine: method II)
- Isooctane/chloroform (honey)
- Hexane (tissue/feed/egg: alternative method, shrimps and fish)

Additional equipment

- Vortex (serum and plasma)
- Bench-centrifuge (serum, plasma, tissue, feed, egg, honey, shrimps, fish)
- Centrifuge refrigerated (milk)
- Evaporation apparatus (serum, plasma, urine: method II, tissue, feed, egg, honey, milk, shrimps, fish)
- pHmeter (urine and water)
- Incubator (urine: method II)
- Balance (tissue, feed, egg, shrimps, fish)
- Homogenizer (tissue, feed, egg, shrimps, fish)
- Rotatory shaker (tissue, feed, egg, honey, shrimps, fish)
- Water bath (tissue, feed, egg, honey, shrimps, fish)
- Micropipette 20-200 μ l with suitable tips
- Multichannel micropipette 50-200 μ l with suitable tips
- ELISA plate or strip reader with a 450nm filter

Storage conditions: +2/+8 °C

Shelf-life: 6 months after manufacture.

Manufacturer

Tecna Srl, ISO9001/UNI EN ISO 9001 – Ed. 2000 certified (SGS, N°. IT01/0291).